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System overview: See Table 1.

Indication

- To prepare the implant socket for common dental implants. If a template is used precaution should be taken so that the instrument does not jam to avoid possible breakage.
- To determine the position, direction and depth of the dental implant.

Clinical sequence

- Determine the correct size of the twist drill and the drilling depth required.
- Determine the size and position of the implant with the aid of radiographs.
- To establish the exact position and the depth of the socket, additional CT scan is recommended.
- In order to prevent any danger to adjacent structures, inspect carefully the area surrounding the preparation site.
- Administer a local anesthesia to the operation site.
- Cut the mucosa and periost to the bone and expose the operation site.
- Mark the position with the a round bur, Apply the pilot bur at the location.
- Drilling should be in a discontinuous motion with constant cooling using saline solution to avoid thermal damage to the bone tissue and to guarantee optimal chip removal.
- Do not allow the drill get jammed or stuck during drilling (risk of fracture). Apply low pressure and drill to the required depth at a speed of 800-1000 rpm. Do not exceed 30 Ncm of torque.
- The progressively larger drills are used to widen the implant site for the implant size selected.

Intended user

The device is allowed to be used by professional dentist.

Patient subject

Healthy adults without serious physical or mental issues.

Contraindication

- Patients with poor bone quality
- Patients with non-controlled periodontal disease or health condition not applicable to dental surgery (such as diabetes or cardiopathy)
- Pregnant women or children are not allowed to implement the surgery.

Side effect

- Perforation of the sinus or mandibular base
- Infections, evident abscesses, fistulas, or inflammation
- Persistent pain, numbness or paresthesia.
- Bone Dehiscence and bone grafting is necessary

WARNING

- The drill is provided in non-sterile condition; therefore it must be cleaned, and sterilized before first time use.
- After use, please immediately proceed to clean and sterilize. Do not let any blood or tissue to dry on the surface after use.
- Immerse the drills in disinfecting solution immediately after use.
- Follow carefully manufacturer's instructions for the use of the cleaning products to prevent any deterioration.
- For recommendations for use (immersion time, concentration, suitability) of the disinfecting and cleaning agents please refer to the manufacturers' instructions.
- Only use non-metal brushes for cleaning to avoid scratches. Do not collide and bend it to avoid scratch and crack.
- Avoid contact with other drills or instruments when cleaning in the ultrasonic bath.

Safety and liability

Cutting hard materials may cause premature blunting; therefore, the drills need to be examined after each use. Replace when necessary.

Precaution:

- Notification of injury caused by the sharp blades of the drill.
- Prevent fall-off and artificial collide and bending for the breaking risk of zirconia dental drill.
- If there is any serious incident occurred in relation to the device, please immediately report to the manufacturer and the competent authority of the member state.
- Safe disposal of the device please refer to Safe Disposal Principal and guideline of medical device.
- Do not use products that exceed the product expiry date.
- Rinse off thoroughly with water and dry carefully.
- Use an autoclave to sterilize these drills following standard procedure regulated by the device's manufacturer.

Cooperated instrument:

Dental drill cooperates certified low-speed dental handpiece.

Prior to operation, the operator undertakes the responsibility to confirm whether the product is suited for the intended uses. The user is responsible for the application of the instrument.

Coho Technology declines any liability for possible damage caused in connection with the use of this instrument.

Sterilization method:

Sterilization is carried out in the autoclave (121-134°C) over a period of 20 minutes as a single cycle. Cleaned and sterilized instruments must be inspected before use, damaged or blunt instruments must be retired. Discard any instruments that have been subjected to any abnormal force (for example dropping to the floor) because the risk of breakage is increased. The operator of these instruments is obligated to make sure that the correct procedure is performed by qualified personnel with appropriate materials and equipment. Work instructions with regard to proper reprocessing of instruments according to EN ISO 17664 Sterilization of medical device.

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