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Zibone Zirconia Dental Implant

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pressure, situations such as peri-implantitis or implant fracture might happen. Patient’s physical condition is critical to successful implant surgery. Therefore, below condition must be considered before implantation:

- Poor bone quality
- Poor oral hygiene
- Health risk due to blood disease or hormone non-homeostasis
- Alcohol or drug abuse
- Implants must be prevented from heavy force after surgery
- Laser engraving is prohibited to implants.
- Therapists are obligated to inform relevant information to patients.

**Precautions：** Surgery for Zibone implants must be executed by professionally trained dentist. Due to highly technical and complex surgical procedure, full surgical evaluation is necessary. Imperfect surgery might lead to functional failure of implant, supporting bone structure loss or other side-effects. Factors such as patients’ health state for implant surgery, patients’ age, or hygiene habits affecting to implants must be fully evaluated.

If there is serious event relevant to the product, please notify immediately to the manufacturer. The device’s safe proposal please refer to safe disposal guidance and proposal of medical device.

### Possible adverse events：

The following adverse events might happen after implant surgery：

- No osseointegration
- Loosening after osseointegration
- Bone Dehiscence and bone grafting is necessary
- Perforation of the sinus or mandibular base
- Infections, evident abscesses, fistulas, or inflammation
- Persistent pain, numbness or paresthesia.
- Excessive bone loss might need surgical intervention.
- Fracture of implant.

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**Sterilization and handling：** The product has been processed via moist heat sterilization. Implants must be stored in original packaging, under clean and dry conditions. If moisture is present in inner package, package damaged, or unintentionally opened, implants must not be used. Before use it, specification number must be checked against package’s label. After taking it out, the product must be prevented to contact with any non-sterilized products. Re-sterilization must be prohibited.

**Surgical Procedure：** Zibone implants must be operated with sterilized instruments in dental operating room. During surgery, drill hole in predetermined location. A minimum space of 1.5mm/2mm should be maintained between implant and implant/natural teeth; a minimum space of 1mm should be maintained between implant and bone edge of lingual side, and buccal side. Preparation of appropriate size and depth for implant site, using torque wrench to thread implant into predetermined depth. Torque force should be within 30-40 Ncm. Implants should be prevented by direct force at initial implantation. Complete osseointegration after eight to twelve weeks, manufacture zirconia prosthetic crown as per general procedure. Surgical drills and product size comparison sheet, please refer to table.

### Notification：

- Please do not throw or press it. Please do not use it after throwing it.
- Because it is a one-piece designed implant, incomplete osseointegration might cause failure easily. Therefore, patients should be tracked for longer period before placing on prosthesis.
- Post-surgery period until osseointegration, overweight occlusion and direct force must be avoided.

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
NOTE: Higher bone density (D1, D2) is suggested going with bone tap.

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